

Product datasheet for **GA208291**

Mouse Pin4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pin4 CRISPRa kit - CRISPR gene activation of mouse protein (peptidyl-prolyl cis/trans isomerase) NIMA-interacting, 4 (parvulin)
Symbol:	Pin4
UniProt ID:	<u>Q9CWW6</u>
Synonyms:	2410002I22Rik; EPVH; Par14
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)