

## Product datasheet for **GA208300**

### Mouse Ttl activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Ttl CRISPRa kit - CRISPR gene activation of mouse tubulin tyrosine ligase
Symbol:	Ttl
UniProt ID:	<a href="#">P38585</a>
Synonyms:	2410003M22Rik; 2700049H19Rik; AI848570
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)