

Product datasheet for **GA208304**

Mouse Pold4 activation kit by CRISPRa

Product data:

| | |
|---------------|----------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits |
| Description: | Pold4 CRISPRa kit - CRISPR gene activation of mouse polymerase (DNA-directed), delta 4 |
| Symbol: | Pold4 |
| UniProt ID: | Q9CWP8 |
| Synonyms: | 2410012M21Rik; A1463381; AW060307; p12; Polds |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)