

Product datasheet for **GA208314**

Mouse Fam174c activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	1600002K03Rik CRISPRa kit - CRISPR gene activation of mouse RIKEN cDNA 1600002K03 gene
Symbol:	Fam174c
Components:	GA208314G1, Fam174c gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208314G2, Fam174c gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208314G3, Fam174c gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9DAZ5</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)