

Product datasheet for **GA208352**

Mouse Ppp1r35 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ppp1r35 CRISPRa kit - CRISPR gene activation of mouse protein phosphatase 1, regulatory subunit 35
Symbol:	Ppp1r35
Components:	GA208352G1, Ppp1r35 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208352G2, Ppp1r35 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208352G3, Ppp1r35 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9D8C8
Synonyms:	2010007H12Rik; 2010011D20Rik; AL024004
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)