

Product datasheet for **GA208358**

Mouse Fam241b activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Fam241b CRISPRa kit - CRISPR gene activation of mouse family with sequence similarity 241, member B
Symbol:	Fam241b
Components:	GA208358G1, Fam241b gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208358G2, Fam241b gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208358G3, Fam241b gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9D882</u>
Synonyms:	AA545165
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)