

Product datasheet for **GA208365**

Mouse Nup43 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Nup43 CRISPRa kit - CRISPR gene activation of mouse nucleoporin 43
Symbol:	Nup43
Components:	GA208365G1, Nup43 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208365G2, Nup43 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208365G3, Nup43 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	P59235
Synonyms:	2610016K01Rik; 2610529I12Rik; AA409950; p42
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)