

Product datasheet for **GA208371**

Mouse Cenps activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Cenps CRISPRa kit - CRISPR gene activation of mouse centromere protein S
Symbol:	Cenps
Components:	GA208371G1, Cenps gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208371G2, Cenps gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208371G3, Cenps gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9D084</u>
Synonyms:	2610040C18Rik; 2810407L01Rik; CENP-S
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)