

## Product datasheet for **GA208400**

### Mouse Dop1b activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Dop1b CRISPRa kit - CRISPR gene activation of mouse DOP1 leucine zipper like protein B
Symbol:	Dop1b
UniProt ID:	<u><a href="#">Q3UHQ6</a></u>
Synonyms:	0610038M01Rik; 2610510B01Rik; AI503246; C21orf5; mKIAA0933
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)