


Product datasheet for **GA208415**

Mouse IgSF23 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	IgSF23 CRISPRa kit - CRISPR gene activation of mouse immunoglobulin superfamily, member 23
Symbol:	IgSF23
Components:	GA208415G1, IgSF23 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208415G2, IgSF23 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208415G3, IgSF23 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	2210010C17Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)