

Product datasheet for **GA208424**

Mouse 2310015D24Rik activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	2310015D24Rik CRISPRa kit - CRISPR gene activation of mouse RIKEN cDNA 2310015D24 gene
Symbol:	2310015D24Rik
Components:	GA208424G1, 2310015D24Rik gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208424G2, 2310015D24Rik gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208424G3, 2310015D24Rik gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)