

Product datasheet for **GA208431**

Mouse Dpf3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Dpf3 CRISPRa kit - CRISPR gene activation of mouse D4, zinc and double PHD fingers, family 3
Symbol:	Dpf3
UniProt ID:	P58269
Synonyms:	2810403B03Rik; 6530402L11Rik; BAF45C; C78788; cer-d4; CERD4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)