


Product datasheet for **GA208433**

Mouse Misp3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	2210011C24Rik CRISPRa kit - CRISPR gene activation of mouse RIKEN cDNA 2210011C24 gene
Symbol:	Misp3
Components:	GA208433G1, Misp3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208433G2, Misp3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208433G3, Misp3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)