

Product datasheet for **GA208447**

Mouse Tmem143 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | Tmem143 CRISPRa kit - CRISPR gene activation of mouse transmembrane protein 143 |
| Symbol: | Tmem143 |
| UniProt ID: | Q8VD26 |
| Synonyms: | 2310076O21Rik; AI481604 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)