

Product datasheet for **GA208451**

Mouse 3300002P13Rik activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	3300002P13Rik CRISPRa kit - CRISPR gene activation of mouse RIKEN cDNA 3300002P13 gene
Symbol:	3300002P13Rik
Components:	GA208451G1, 3300002P13Rik gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208451G2, 3300002P13Rik gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208451G3, 3300002P13Rik gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)