

Product datasheet for **GA208463**

Mouse 2310068J16Rik activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	2310068J16Rik CRISPRa kit - CRISPR gene activation of mouse RIKEN cDNA 2310068J16 gene
Symbol:	2310068J16Rik
Components:	GA208463G1, 2310068J16Rik gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208463G2, 2310068J16Rik gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208463G3, 2310068J16Rik gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)