

Product datasheet for **GA208478**

Mouse Cd3eap activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Cd3eap CRISPRa kit - CRISPR gene activation of mouse CD3E antigen, epsilon polypeptide associated protein
Symbol:	Cd3eap
Components:	GA208478G1, Cd3eap gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208478G2, Cd3eap gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208478G3, Cd3eap gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q76KJ5
Synonyms:	1190028F09; 2610103M17Rik; ASE-1; Ase1; CAST; PAF49
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)