

## Product datasheet for **GA208483**

### Mouse Ppp4r1 activation kit by CRISPRa

#### Product data:

|               |  |
|---------------|--|
| Product Type: | CRISPRa Kits   |
| Description:  | Ppp4r1 CRISPRa kit - CRISPR gene activation of mouse protein phosphatase 4, regulatory subunit 1 |
| Symbol:       | Ppp4r1   |
| UniProt ID:   | <a href="#">Q8K2V1</a>   |
| Synonyms:     | 3110001J10Rik; Pp4r1   |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                          |

#### Product images:



[View online >](#)