


Product datasheet for **GA208500**

Mouse Kctd2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Kctd2 CRISPRa kit - CRISPR gene activation of mouse potassium channel tetramerisation domain containing 2
Symbol:	Kctd2
Components:	GA208500G1, Kctd2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208500G2, Kctd2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208500G3, Kctd2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	2310012I15Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)