

## Product datasheet for **GA208515**

### Mouse Tspan15 activation kit by CRISPRa

#### Product data:

|               |   |
|---------------|---|
| Product Type: | CRISPRa Kits  |
| Description:  | Tspan15 CRISPRa kit - CRISPR gene activation of mouse tetraspanin 15    |
| Symbol:       | Tspan15   |
| UniProt ID:   | <a href="#">F7BWT7</a>  |
| Synonyms:     | 1110036D12Rik; 2700063A19Rik; AW048364; Tm4sf15                         |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

#### Product images:



[View online >](#)