

Product datasheet for **GA208523**

Mouse **Taf15** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	Taf15 CRISPRa kit - CRISPR gene activation of mouse TATA-box binding protein associated factor 15
Symbol:	Taf15
Components:	GA208523G1, Taf15 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208523G2, Taf15 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208523G3, Taf15 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	68kDa; 2610111C21Rik; Taf2n; TAFII68
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)