

Product datasheet for **GA208524**

Mouse Mir124-2hg activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Mir124-2hg CRISPRa kit - CRISPR gene activation of mouse Mir124-2 host gene (non-protein coding)
Symbol:	Mir124-2hg
Components:	GA208524G1, Mir124-2hg gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208524G2, Mir124-2hg gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208524G3, Mir124-2hg gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)