

## Product datasheet for **GA208525**

### Mouse Cd248 activation kit by CRISPRa

#### Product data:

|               |   |
|---------------|---|
| Product Type: | CRISPRa Kits  |
| Description:  | Cd248 CRISPRa kit - CRISPR gene activation of mouse CD248 antigen, endosialin |
| Symbol:       | Cd248   |
| UniProt ID:   | <a href="#">Q91V98</a>  |
| Synonyms:     | 2610111G01Rik; AI842296; Cd164I1; Tem1  |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)       |

#### Product images:



[View online >](#)