

Product datasheet for **GA208533**

Mouse Ckap2l activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ckap2l CRISPRa kit - CRISPR gene activation of mouse cytoskeleton associated protein 2-like
Symbol:	Ckap2l
Components:	GA208533G1, Ckap2l gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208533G2, Ckap2l gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208533G3, Ckap2l gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q7TS74</u>
Synonyms:	2010016H04Rik; 2610318C08Rik; AV070319; radmis; rdm
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)