


Product datasheet for **GA208545**

Mouse Rnf167 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Rnf167 CRISPRa kit - CRISPR gene activation of mouse ring finger protein 167
Symbol:	Rnf167
Components:	GA208545G1, Rnf167 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208545G2, Rnf167 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208545G3, Rnf167 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q91XF4</u>
Synonyms:	0610010G05Rik; 5730408C10Rik; AU019305; AV328608
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)