

Product datasheet for **GA208555**

Mouse Tmtc4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmtc4 CRISPRa kit - CRISPR gene activation of mouse transmembrane and tetratricopeptide repeat containing 4
Symbol:	Tmtc4
Components:	GA208555G1, Tmtc4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208555G2, Tmtc4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208555G3, Tmtc4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8BG19
Synonyms:	4930403J22Rik; 5730419O14Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)