

Product datasheet for **GA208568**

Mouse 5730455P16Rik activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	5730455P16Rik CRISPRa kit - CRISPR gene activation of mouse RIKEN cDNA 5730455P16 gene
Symbol:	5730455P16Rik
Components:	GA208568G1, 5730455P16Rik gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208568G2, 5730455P16Rik gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208568G3, 5730455P16Rik gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9CYI0</u>
Synonyms:	RP23-252E4.2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)