

Product datasheet for **GA208591**

Mouse *Sik3* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Sik3 CRISPRa kit - CRISPR gene activation of mouse SIK family kinase 3
Symbol:	Sik3
UniProt ID:	Q6P4S6
Synonyms:	5730525O22Rik; 9030204A07Rik; AI447252; AI846254; AI849445; mKIAA0999; Qsk; SIK-3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)