

Product datasheet for **GA208600**

Mouse **Gpr137c** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	Gpr137c CRISPRa kit - CRISPR gene activation of mouse G protein-coupled receptor 137C
Symbol:	Gpr137c
Components:	GA208600G1, Gpr137c gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208600G2, Gpr137c gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208600G3, Gpr137c gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	E9Q343
Synonyms:	6330416L11Rik; Gm908; TM7SF1L2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)