

## Product datasheet for GA208601

### Mouse Medag activation kit by CRISPRa

#### Product data:

|                      |                                                                                                                                                                                                                                                                                                            |
|----------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type:</b> | CRISPRa Kits                                                                                                                                                                                                                                                                                               |
| <b>Description:</b>  | Medag CRISPRa kit - CRISPR gene activation of mouse mesenteric estrogen dependent adipogenesis                                                                                                                                                                                                             |
| <b>Symbol:</b>       | Medag                                                                                                                                                                                                                                                                                                      |
| <b>Components:</b>   | GA208601G1, Medag gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA208601G2, Medag gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA208601G3, Medag gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| <b>UniProt ID:</b>   | <u><a href="#">Q14BA6</a></u>                                                                                                                                                                                                                                                                              |
| <b>Synonyms:</b>     | 6330406I15Rik; Awms3; mAWMS3; MEDA-4; Meda4                                                                                                                                                                                                                                                                |
| <b>Format:</b>       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                    |

#### Product images:

