

## Product datasheet for **GA208622**

### Mouse Ankib1 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                               |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                  |
| Description:  | Ankib1 CRISPRa kit - CRISPR gene activation of mouse ankyrin repeat and IBR domain containing 1                                                                                                                                                                                                               |
| Symbol:       | Ankib1                                                                                                                                                                                                                                                                                                        |
| Components:   | GA208622G1, Ankib1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA208622G2, Ankib1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA208622G3, Ankib1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u>Q6ZPS6</u>                                                                                                                                                                                                                                                                                                 |
| Synonyms:     | 2310061P20Rik; 4631416I11Rik; AW494740; C80642; mKIAA1386                                                                                                                                                                                                                                                     |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                       |

#### Product images:



[View online »](#)