


Product datasheet for **GA208636**

Mouse P2ry12 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	P2ry12 CRISPRa kit - CRISPR gene activation of mouse purinergic receptor P2Y, G-protein coupled 12
Symbol:	P2ry12
UniProt ID:	Q9CPV9
Synonyms:	2900079B22Rik; 4921504D23Rik; P2Y12
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)