


Product datasheet for **GA208667**

Mouse Cdkn2aip activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Cdkn2aip CRISPRa kit - CRISPR gene activation of mouse CDKN2A interacting protein
Symbol:	Cdkn2aip
Components:	GA208667G1, Cdkn2aip gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208667G2, Cdkn2aip gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208667G3, Cdkn2aip gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8BI72
Synonyms:	492151116Rik; AW208986; CARF
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)