

Product datasheet for **GA208673**

Mouse Mmrn1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Mmrn1 CRISPRa kit - CRISPR gene activation of mouse multimerin 1
Symbol:	Mmrn1
Components:	GA208673G1, Mmrn1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208673G2, Mmrn1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208673G3, Mmrn1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>B2RPV6</u>
Synonyms:	4921530G03Rik; Emilin4; Mmrn
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)