

Product datasheet for **GA208693**

Mouse Spef1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Spef1 CRISPRa kit - CRISPR gene activation of mouse sperm flagellar 1
Symbol:	Spef1
Components:	GA208693G1, Spef1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208693G2, Spef1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208693G3, Spef1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q99JL1
Synonyms:	4931426K16Rik; Clamp
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)