

## Product datasheet for **GA208728**

### Mouse Tssk4 activation kit by CRISPRa

#### Product data:

|               |   |
|---------------|---|
| Product Type: | CRISPRa Kits  |
| Description:  | Tssk4 CRISPRa kit - CRISPR gene activation of mouse testis-specific serine kinase 4 |
| Symbol:       | Tssk4   |
| UniProt ID:   | <a href="#">Q9D411</a>  |
| Synonyms:     | 1700020B19Rik; 4933424F08Rik; TSK-4; TSSK-4   |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)             |

#### Product images:



[View online >](#)