


Product datasheet for **GA208739**

Mouse Tmem217 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmem217 CRISPRa kit - CRISPR gene activation of mouse transmembrane protein 217
Symbol:	Tmem217
Components:	GA208739G1, Tmem217 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208739G2, Tmem217 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208739G3, Tmem217 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	4933413N12Rik; EG622644
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)