

Product datasheet for **GA208745**

Mouse Eri2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Eri2 CRISPRa kit - CRISPR gene activation of mouse exoribonuclease 2
Symbol:	Eri2
Components:	GA208745G1, Eri2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208745G2, Eri2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208745G3, Eri2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q5BKS4
Synonyms:	4933424N09Rik; Exod1; mKIAA1504
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)