


Product datasheet for **GA208768**

Mouse Daw1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Daw1 CRISPRa kit - CRISPR gene activation of mouse dynein assembly factor with WDR repeat domains 1
Symbol:	Daw1
Components:	GA208768G1, Daw1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208768G2, Daw1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208768G3, Daw1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	4930563E19Rik; 4933429D11Rik; b2b1116Clo; b2b1584Clo; Wdr69
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)