

Product datasheet for **GA208793**

Mouse Mfap3l activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Mfap3l CRISPRa kit - CRISPR gene activation of mouse microfibrillar-associated protein 3-like
Symbol:	Mfap3l
Components:	GA208793G1, Mfap3l gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208793G2, Mfap3l gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208793G3, Mfap3l gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9D3X9
Synonyms:	4933428A15Rik; 5430405D20Rik; AI461995; AW125052; mKIAA0626; NYD-sp9
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)