

Product datasheet for **GA208806**

Mouse Wdr31 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Wdr31 CRISPRa kit - CRISPR gene activation of mouse WD repeat domain 31
Symbol:	Wdr31
Components:	GA208806G1, Wdr31 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208806G2, Wdr31 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208806G3, Wdr31 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9JHB4
Synonyms:	5430402I10Rik; spWD
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)