

Product datasheet for **GA208811**

Mouse Pdss2 activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | Pdss2 CRISPRa kit - CRISPR gene activation of mouse prenyl (solanesyl) diphosphate synthase, subunit 2 |
| Symbol: | Pdss2 |
| UniProt ID: | Q33DR3 |
| Synonyms: | 5430420P03Rik; Gm60; kd; mDLP1; Plmp |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)