


## Product datasheet for **GA208819**

### Mouse Chd6 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                         |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                            |
| Description:  | Chd6 CRISPRa kit - CRISPR gene activation of mouse chromodomain helicase DNA binding protein 6                                                                                                                                                                                                          |
| Symbol:       | Chd6                                                                                                                                                                                                                                                                                                    |
| Components:   | GA208819G1, Chd6 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA208819G2, Chd6 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA208819G3, Chd6 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <a href="#">A3KFM7</a>                                                                                                                                                                                                                                                                                  |
| Synonyms:     | 5430439G14Rik; 6330406J24Rik; CHD-6                                                                                                                                                                                                                                                                     |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                 |

#### Product images:



[View online »](#)