

## Product datasheet for **GA208823**

### Mouse Cdc3 activation kit by CRISPRa

#### Product data:

<b>Product Type:</b>	CRISPRa Kits
<b>Description:</b>	5430419D17Rik CRISPRa kit – CRISPR gene activation of mouse RIKEN cDNA 5430419D17 gene
<b>Symbol:</b>	Cdc3
<b>Components:</b>	GA208823G1, Cdc3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208823G2, Cdc3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208823G3, Cdc3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
<b>UniProt ID:</b>	<u><a href="#">Q8BZE1</a></u>
<b>Format:</b>	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:

