

Product datasheet for **GA208847**

Mouse Vps35l activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Vps35l CRISPRa kit - CRISPR gene activation of mouse VPS35 endosomal protein sorting factor like
Symbol:	Vps35l
Components:	GA208847G1, Vps35l gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208847G2, Vps35l gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208847G3, Vps35l gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8BWQ6
Synonyms:	2610029J22Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)