

Product datasheet for **GA208864**

Mouse Pogk activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pogk CRISPRa kit - CRISPR gene activation of mouse pogo transposable element with KRAB domain
Symbol:	Pogk
UniProt ID:	<u>Q80TC5</u>
Synonyms:	9130401E23Rik; BASS2; mKIAA1513
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)