

Product datasheet for **GA208888**

Mouse Sanbr activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	0610010F05Rik CRISPRa kit – CRISPR gene activation of mouse RIKEN cDNA 0610010F05 gene
Symbol:	Sanbr
Components:	GA208888G1, Sanbr gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208888G2, Sanbr gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208888G3, Sanbr gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q68FF0</u>
Synonyms:	Kiaa1841; mKIAA1841
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

