

Product datasheet for **GA208907**

Mouse Dhx35 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Dhx35 CRISPRa kit - CRISPR gene activation of mouse DEAH (Asp-Glu-Ala-His) box polypeptide 35
Symbol:	Dhx35
Components:	GA208907G1, Dhx35 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208907G2, Dhx35 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208907G3, Dhx35 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	1200009D07Rik; Ddx35
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)