

## Product datasheet for **GA208912**

### Mouse Dhx34 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                            |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                               |
| Description:  | Dhx34 CRISPRa kit - CRISPR gene activation of mouse DEAH (Asp-Glu-Ala-His) box polypeptide 34                                                                                                                                                                                                              |
| Symbol:       | Dhx34                                                                                                                                                                                                                                                                                                      |
| Components:   | GA208912G1, Dhx34 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA208912G2, Dhx34 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA208912G3, Dhx34 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q9DBV3</a></u>                                                                                                                                                                                                                                                                              |
| Synonyms:     | 1200013B07Rik; 1810012L18Rik; Ddx34; mKIAA0134                                                                                                                                                                                                                                                             |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                    |

#### Product images:



[View online »](#)