

Product datasheet for **GA208914**

Mouse Stk11ip activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Stk11ip CRISPRa kit - CRISPR gene activation of mouse serine/threonine kinase 11 interacting protein
Symbol:	Stk11ip
Components:	GA208914G1, Stk11ip gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208914G2, Stk11ip gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208914G3, Stk11ip gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q3TAA7</u>
Synonyms:	1200014D22Rik; BB131189; LIP1; LKB1IP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)